



**Note:** Please contact CSAN® at 1 (800) 267-2726 if your device issue is not solved by the User Guide.

**1**

**DEVICE CONNECTED TO POWER**

When you plug the CSAN® Pronto® AC adapter into the power source, the LED light will turn yellow



**2**

**DEVICE CONNECTED TO NETWORK**

When you connect the CSAN® Pronto® to the internet, the LED light will turn green





Note: If the LED light is green but the device is unavailable to run a test, contact CSAN® as it is likely that the network is blocked (firewall).



**3**

**DEVICE READY TO RUN A TEST**

When you insert the test strip into the CSAN® Pronto®, the LED light will turn blue. You can then click Begin Test within the CSAN Patient Care Portal®





**Note:** Please contact CSAN® at 1 (800) 267-2726 if your device issue is not solved by the User Guide.

## 4 DEVICE RUNNING A TEST

When the CSAN® Pronto® is running the test, the LED light will turn white.



## 5 DEVICE ISSUE

If the CSAN® Pronto® LED light turns red there is likely a connectivity issue or an issue with the device itself. Verify the internet connection. If this does not resolve, contact CSAN® at 1-800-267-2726 for assistance.

