



Note: Please contact CSAN® at 1 (800) 267-2726 if your device problem is not solved by the User Guide.

1 HOW TO TURN YOUR IOS DEVICES INTO A WIFI HOTSPOT

To set up a personal hotspot on your iPhone or iPad (WiFi + Cellular), go to **Settings > Personal Hotspot** and toggle it on. Your device will broadcast its own wireless network and name.

Take note of the WiFi password.

Download the CSAN Pronto APP



IOS 8.0 or later versions




2 WIFI HOTSPOT WITH ANDROID DEVICES

You can either pull down and tap the Mobile Hotspot option in the menu to activate. Or Navigate to **Connections > Mobile Hotspot and Tethering > Mobile Hotspot** and toggle it on. A pop-up window will appear, where you can find the WiFi password for devices that want to connect.

Download the CSAN® Pronto® APP



Android 5.1 and up



Tip! The steps to connect to a WiFi network vary depending on your device and software. See your device's User Guide for more information.

3 HOW TO CONNECT YOUR CSAN® PRONTO® DEVICE TO YOUR WIFI HOTSPOT




Access the CSAN® Pronto® App

Connect to your device via Bluetooth

Select the desired device

Scan for wireless networks

From a different device from the hotspot





4 WiFi HOTSPOT CONNECTION


Select your Personal Hotspot




Enter your WiFi Hotspot password




Enter Credentials



Connecting



CSAN[®] Pronto[®] App should be Connected



The LED light of the device will turn GREEN



- **Hotspot connection** is the same as WiFi connection, except the internet is provided by a cellular company instead of an internet provider in your home or office.
- To ensure proper functioning, make sure your hotspot device has a **strong signal** and **LTE**. Note that cellular signal strength may vary from one area to the other.
- The **CSAN[®] Pronto[®] App** is a one-time requirement to make the connection to WiFi, you will not require this App to run tests. Its only purpose is to establish a connection for your CSAN[®] Pronto[®] device. **If multiple users** are using the same CSAN[®] Pronto[®] device, you may need to reestablish a hotspot connection with the App.
- **Once you have established a connection**, access the CSAN Patient Care Portal[®] : demo portal at csan.ca/demo or the production portal at www.clozaril.ca.

Before leaving the clinic with your mobile device:

- Turn on your mobile device hotspot connection and **leave the screen open** on that window.
- Use a **separate device** (i.e. tablet or mobile phone) to go on the CSAN[®] Pronto[®] APP and scan for devices. Once the list appears, select the desired device.
- Scan for wireless networks and **select your personal hotspot** WiFi connection.
- Once connected to the CSAN[®] Pronto[®] device (GREEN), you can turn off the other device used for connection and leave the clinic with your hotspot connected mobile device to execute your CSAN[®] Pronto[®] tests throughout the day.

At the blood collection site:

- Turn on your **personal hotspot** on your mobile device and remain on this screen.
- Plug the CSAN[®] Pronto[®] device into a **power source** and wait for the **LED light** of the device to **turn GREEN**.
- Open the **CSAN Patient Care Portal[®]** via the internet browser of your mobile device and sign in with your credentials.
- Proceed with the usual steps to **execute a CSAN[®] Pronto[®] test**.
- Once the test is complete, we recommend to turn-off the mobile hotspot on your device.

Tips!

- **Once back from the road**, to ensure proper software update, please reconnect your CSAN[®] Pronto[®] device to your standard WiFi network and to an electrical outlet.
- **If you are sharing the usage of your CSAN[®] Pronto[®] device with colleagues**, to limit future Bluetooth pairing issues, unpair your mobile device with your CSAN[®] Pronto[®] device.
- **Recharge** all of your complimentary mobile devices for future use.